

Shining[®] Scan Checklist for Creating Orders

Fill out the **Patient Name** (create a new record if necessary)

Select the correct **Dentistry Type** to determine workflow

- Restorative - for crown preps
 - Implants - for scan bodies
 - Removable Denture - for denture type cases
 - Orthodontic - for clear aligners
 - Other - for extra scans (ex. surgical guides)
-

Select the correct tooth sites and configure accordingly

- Correct material
 - For restorative orders, set Pre-op to Yes
 - Select a shade (if using the Vita Shade guide) or enter shade into the notes
-

*** Enter a return date in the **Notes field** ***