

## Shining® Scan Checklist for Creating Orders

---

☐

Fill out the **Patient Name** (create a new record if necessary)

.....

☐

Select the correct **Dentistry Type** to determine workflow

- Restorative - for crown preps
  - Implants - for scan bodies
  - Removable Denture - for denture type cases
  - Orthodontic - for clear aligners
  - Other - for extra scans (ex. surgical guides)
- .....

☐

Select the correct tooth sites and configure accordingly

- Correct material
  - For restorative and implant orders, check the box for a Pre-op
  - Select a shade (if using the Vita Shade guide) or enter shade into the notes
- .....

☐

\*\*\* Enter a return date in the **Notes field** \*\*\*